

DESCRIPTION OF PRODUCT: A Pelvic Holder to prevent sliding. For chair application only.

Rx ONLY



INDICATIONS FOR USE:

- Patients assessed to be at risk of injury from a fall.
- Patients needing a positioning device for added safety while in a chair.
- Patients who have a tendency to slide down in a chair.

CONTRAINDICATIONS:

- **DO NOT** use on a patient who is or becomes highly aggressive, combative, agitated, or suicidal.
- **DO NOT** use on patients with: ostomy, colostomy, or G-tubes; hernias, severe Cardio Obstructive Pulmonary Disease (COPD); or with post-surgery tubes, incisions, catheters, or monitoring lines. These could be disrupted by a restraint.

ADVERSE REACTIONS

- Severe emotional, psychological, or physical problems may occur: if the applied device is uncomfortable; or if it severely limits movement. If symptoms of these problems ever appear for any reason, get help from a qualified medical authority and find a less restrictive, product or intervention.



POSEY PELVIC HOLDERS

- REF 4427** Breezeline Mesh, with buckles
- REF 4430** Padded with foam liner and quick-release buckles

APPLICATION INSTRUCTIONS:

⚠ WARNING Make sure patient wears proper undergarments to protect skin.

⚠ CAUTION Before use, check device for damage. Discard if you have any questions about patient safety.

1. Lay the pelvic holder on the chair with the narrow side to the back of the chair. The label will be face down in the seat (fig. 1).




2. Bring the ends of the connecting straps on the narrow end, down between the seat and the wheelchair sides.

3. Secure the straps behind the backrest.

4. Position the patient as far back in the seat as possible, with the buttocks against the back of the chair.

5. Bring the wide part of the pelvic holder up between the patient's legs.

6. Bring the ends of the connecting straps on the wide end, down between the seat and wheelchair sides at a 45-degree angle (fig. 2).



7. Secure the straps behind the backrest as low as possible (fig. 3).

8. Check that the straps are secure and will not change position, loosen, or tighten if the patient pulls on them, or if the chair is tilted or adjusted.



⚠ WARNING

Heed these warnings to reduce the risk of serious injury or death:

- Monitor skin conditions in the groin area frequently. If the patient slides down or forward, pelvic straps may damage the skin.
- There is a risk of chest compression or suffocation, if the patient's body weight is suspended off the chair seat. Use extreme caution with chair cushions. If a cushion dislodges, straps may loosen and allow the patient to slide off the seat (fig. 4).
- Monitor per facility policy to ensure that the patient cannot slide down, or fall off the chair seat and become suspended (fig. 4).



FIG. 4

- **STOP USE AT ONCE:** if the patient has a tendency to slide forward or down in the device; or is able to self-release.

ADDITIONAL SAFETY AND LAUNDERING INSTRUCTIONS ON OTHER SIDE



