

# Torso Support 3656 CHECKLIST



Facility		Contact	
Facility Name		Name	
Street Address		Title	
City/State/Zip		Nursing Unit	
Instructor		Date	
Signature			

Description	Competency Level	Date	Comments
Competency Level Definition - (1) Needs remediation and re-evaluation, (2) Performs skill with accuracy under supervision, (3) Performs skill with accuracy independently			
<b>3656S</b> 48"L x 5"W 122 cm x 13 cm (Waist Strap) 34"L x 3"W 86 cm x 8 cm (Shoulder Strap)			
<b>3656M</b> 54"L x 5"W 137 cm x 13 cm (Waist Strap) 37"L x 3"W 94 cm x 8 cm (Shoulder Strap)			
<b>3656L</b> 60"L x 5"W 152 cm x 13 cm (Waist Strap) 41"L x 3"W 104 cm x 8 cm (Shoulder Strap)			
<b>3656XL</b> 72"L x 5"W 183 cm x 13 cm (Waist Strap) 41"L x 3"W 104 cm x 8 cm (Shoulder Strap)			
<b>3656XXL</b> 84"L x 5"W 213 cm x 13 cm (Waist Strap) 46"L x 3"W 117 cm x 8 cm (Shoulder Strap)			
<b>APPLICATION INSTRUCTIONS:</b> 1. Seat the patient with the hips against the chair back.			
2. Wrap the chest strap around the patient and the chair. Press hook and loop firmly together. The patient should be comfortably seated and breathing easily.			

Rx ONLY

Before use, make sure to read the instructions accompanying the product. A copy of the instruction sheet is included with the product and can also be downloaded at [www.posey.com](http://www.posey.com).

5635 Peck Road, Arcadia, CA 91006-0020 USA • Phone: 1.800.447.6739 • Fax: 1.800.767.3933 • [www.posey.com](http://www.posey.com)



Description	Competency Level	Date	Comments
Competency Level Definition - (1) Needs remediation and re-evaluation, (2) Performs skill with accuracy under supervision, (3) Performs skill with accuracy independently			
3. Bring the shoulder straps over the shoulders, crossing them in an "X" behind the seat rest. Secure the shoulder straps to chest strap.			
<b>Note:</b> The "hook and loop" may be secured in front of the patient for "self-release," or behind the patient for "assisted-release". The proper medical authority should determine which way the belt is to be utilized.			