

# Twice-As-Tough™ Cuffs for Stretchers 2750, 2755 CHECKLIST



Facility		Contact	
Facility Name		Name	
Street Address		Title	
City/State/Zip		Nursing Unit	
Instructor		Date	
Signature			

Description	Competency Level	Date	Comments
Competency Level Definition - (1) Needs remediation and re-evaluation, (2) Performs skill with accuracy under supervision, (3) Performs skill with accuracy independently			
<b>2750</b> Wrist pair, regular stretcher (blue neoprene) <b>2755</b> Ankle pair, regular stretcher (red neoprene)			
<b>APPLICATION INSTRUCTIONS:</b> <b>⚠ WARNING</b> <b>Connecting Straps</b> (repeat steps 1-4 for each side) Before each use, check cuffs and straps for cracks, tears, and/or excessive wear or stretch; cracked or broken buckles or locks; and/or that hook and loop adheres securely, as these may allow patient to remove cuff. Discard if device is damaged.  If you have any questions about patient safety, ask the appropriate medical authority for alternatives.			
1. Place the cuffs on the stretcher at the wrist or ankle position.			
2. Bring the ends of the strap down through the inside of the side rails so they do not interfere with the movement of the side rails. Thread the end of the strap over the top, around the frame, and connect the blue or red hook and loop together.			

Rx ONLY

Before use, make sure to read the instructions accompanying the product. A copy of the instruction sheet is included with the product and can also be downloaded at [www.posey.com](http://www.posey.com).

Posey Company • 5635 Peck Road, Arcadia, CA 91006-0020 USA • Phone: 1.800.447.6739 • Fax: 1.800.767.3933 • [www.posey.com](http://www.posey.com)



Description	Competency Level	Date	Comments
Competency Level Definition - (1) Needs remediation and re-evaluation, (2) Performs skill with accuracy under supervision, (3) Performs skill with accuracy independently			
3. Pull the strap snug.			
4. Position the patient on the stretcher.			
<b>Applying the Cuffs</b> (repeat steps 1-2 for each limb) 1. Wrap the neoprene piece (the blue or red side should be positioned against the skin) around the wrist/ankle. Attach the black hook and loop pieces together, followed by the blue or red hook and loop pieces. The fuzzy piece should be sandwiched between the two pieces of hook. Be sure to overlap at least one inch (3 cm).			
2. Press the hook and loop closure together firmly and make sure it adheres securely. Slide ONE finger (flat) between the cuff and the inside of the patient's wrist/ankle to ensure proper fit. The cuffs must be snug enough to prevent escape, but not interfere with circulation.			