

# THE PoseyBED

## Model 8070

The Posey Bed is part of the Fall Management Protocol and is used in more than 600 Hospitals, LTC and Homecare sites worldwide.




- *Indicated for behaviorally, physically or mentally compromised patients at risk of serious injury from a fall.*
- *A less restrictive alternative to chemical and physical protective restraints.*

# THE Posey BED Model 8070

U.S. Patent No. 7,047,991

1. Medical Bed. Fully automatic low hospital bed (13 ½" floor to mattress pan).
2. Patented Enclosed Mattress Compartment. Helps reduce the risk of patient entrapment.
3. Nylon "U" Side Panels with a single zipper on each panel for easy patient access.
4. Head and Foot Access Panels with a single zipper on each panel.
5. Quick-Release Buckles clip to zipper pull-tabs for added patient security.
6. Perimeter Guard. These function like side rails for temporary use during patient care and must be unzipped when side panels are closed for easier patient access and better patient visibility.
7. Side Panel Tube Ports are located on each side at the foot of the bed. Use these to run tubes into the patient compartment.
8. Locking Casters secure bed during patient transfer to and from the bed and when patient is unattended.
9. Transfer Brakes. These engage when bed is in the lowest position. They help add stability and prevent lateral movement during patient transfer.
10. Drainage Ports are located on each side of the head of the bed.



Please detach and fax to 1.626.443.5014, or call us at 1.800.44.POSEY

Name \_\_\_\_\_

Title \_\_\_\_\_

Facility \_\_\_\_\_

Address \_\_\_\_\_

City \_\_\_\_\_ State \_\_\_\_\_ Zip \_\_\_\_\_

Telephone \_\_\_\_\_ Best time to call \_\_\_\_\_ am \_\_\_\_\_ pm

Email \_\_\_\_\_

**yes!**

Please have a Posey Sales Representative call me to discuss the Posey Bed Model 8070.



Fax or mail to:



Posey Company  
 Attn: Marketing Dept.  
 5635 Peck Road • Arcadia, CA 91006-0020 USA  
 Tel: 1.800.44.POSEY (447.6739) • Fax 1.626.443.5014 • www.posey.com