


## TECHNICAL SPECIFICATIONS



### DIMENSIONS:

**Area:** 24 square feet - **Volume (Inside):** 70 cubic feet (2 cubic meters)

**Total Weight completely assembled (bed frame and canopy):** 278 pounds (126 Kg)

**Total Weight of the canopy only:** 17 pounds (8 Kg)

**Maximum Back Angle:** 70° to Horizontal

**Maximum Knee Angle:** 35° to Horizontal

**Maximum Knee to Leg Angle:** 35° to Horizontal

### BED ELECTRICAL SPECIFICATIONS:

- **Power:** 108-132 VAC
- **Frequency:** 60Hz
- **Current Rating:** 1.6 AMP (110-120V Unit)
- **Classification:** Class 2, Type B (Double Insulated)
- **Circuit Protection:** 250V, F1.6A (110-120V Unit)
- **Mode of Operation:** 10% Maximum duty cycle