




Posey® Lateral Wedges

Posey Lateral wedges are lightweight, polystyrene filled wedges that offer both positioning support and comfort. Lateral wedges are ideal for keeping bed-ridden patients repositioned to eliminate pressure points over bony prominences and thereby reducing the incidence of pressure ulcers.

Posey's unique construction of an outer VINYL MESH, breathable fabric shell and a Non-absorbent, POLYSTYRENE BEAD fill allows air to flow through the wedge thereby eliminating discomfort often associated with conventional foam pillows, which do not breath and build up heat. The conformability of the bead fill design offers a better solution than pillows as Posey Lateral Wedges stay better positioned and will not compress like foam.



#6309

Each lateral wedge features a zippered closure to allow the addition or removal of beads, as it may be desirable to change the amount of fill to better accommodate certain clients.

Posey Lateral Wedges are also more cost effective than pillows as they can reduce nursing labor costs. With the lightweight Posey wedge it is possible to turn the patient without assistance. The Posey Lateral Wedge is easily placed in position unlike foam pillows, which compress and are difficult to keep in place.

Laundering Instructions:

Machine washable in warm temperature. Air dry.



Description	Typical Client	Dimensions	Posey Item #	Replaces MEDCO Item #
Lateral Wedge - Small	Pediatrics or for cervical use	14½" L x 4" W x 5½"H	REF #6309S	1LW2-S
Lateral Wedge - Medium	Patients under 180 lbs	16" L x 6" W x 7" H	REF #6309M	1LW3-M
Lateral Wedge - Large	Patients over 180 lbs	19" L x 7" W x 9" H	REF #6309L	1LW4-L
Lateral Wedge - Extra Long	Very Tall Patients	23" L x 7" W x 9" H	REF #6309XL	1LW5-XL
Beads Refill Kit	N/A		REF #8124	3PK

Materials used meet the requirements of FMVSS-302 and CA #117. Section E for Flammability

J.T. Posey Company

5635 Peck Road • Arcadia, CA 91006-0020 USA • Tel: 800-447-6739 or 626-443-3143 • Fax: 800-767-3933 or 626-443-5014 • www.posey.com
 © 2004 J.T. Posey Company. All rights reserved.



MDSS
 Burckhardtstr. 1,
 30163, Hannover, Germany



M6141 123004