



# Posey® Leg & Foot Hugger™

Leg & Foot Hugger provides cushioned leg support and helps prevent resident's feet from slipping off footrest and getting caught under chair. Footrest extender helps to prevent footdrop by increasing the base of support and elevating the footrest. Also helps keep resident's feet from slipping off the footrests. Cushioned with foam padding and covered in water-repellent, wipe-clean vinyl. Fits 18" wide wheelchair. All models meet CA #117 flame retardancy standards.

- Cat. No. 6351**      **Foot Hugger – 1" footrest extender**
- Cat. No. 6352**      **Foot Hugger – 2" footrest extender**
- Cat. No. 6353**      **Foot Hugger – 3" footrest extender**
- Cat. No. 6354**      **Foot Hugger – 4" footrest extender**
- Cat. No. 6355**      **Leg Hugger – no footrest extender**

### Application Instructions:

1. Position the patient in the wheelchair and adjust the footrests to fit.
2. Release the quick-release buckles on the attachment straps.
3. Position the Hugger **between** the footrest bars with the extender resting on top of the footrests as shown in the photo.
4. Wrap the top straps around the outside of each footrest bar and fasten the quick-release buckles. Follow same procedure for the lower straps.
5. Adjust the straps so the leg rest sits **flush** between the footrest bars. This will allow more room on the footrests for the resident's feet.

NOTE: Extenders are available in 1", 2", 3", and 4" thicknesses and may be purchased separately. Call Customer Service at 800-447-6739.



6351

**J.T. Posey Company** 5635 Peck Road • Arcadia, CA 91006-0020 USA • Tel: 800-447-6739 or 626-443-3143 • Fax: 800-767-3933 or 626-443-5014

EU AUTH. PERSON: (MDD 93/42/EEC) MDSS GmbH, Burckhardtstr. 1, D-30163, Hannover, Germany

© 1996 J.T. Posey Company. All rights reserved.



M5015 071598



# Posey® Leg & Foot Hugger™

Leg & Foot Hugger provides cushioned leg support and helps prevent resident's feet from slipping off footrest and getting caught under chair. Footrest extender helps to prevent footdrop by increasing the base of support and elevating the footrest. Also helps keep resident's feet from slipping off the footrests. Cushioned with foam padding and covered in water-repellent, wipe-clean vinyl. Fits 18" wide wheelchair. All models meet CA #117 flame retardancy standards.

- Cat. No. 6351**      **Foot Hugger – 1" footrest extender**
- Cat. No. 6352**      **Foot Hugger – 2" footrest extender**
- Cat. No. 6353**      **Foot Hugger – 3" footrest extender**
- Cat. No. 6354**      **Foot Hugger – 4" footrest extender**
- Cat. No. 6355**      **Leg Hugger – no footrest extender**

### Application Instructions:

1. Position the patient in the wheelchair and adjust the footrests to fit.
2. Release the quick-release buckles on the attachment straps.
3. Position the Hugger **between** the footrest bars with the extender resting on top of the footrests as shown in the photo.
4. Wrap the top straps around the outside of each footrest bar and fasten the quick-release buckles. Follow same procedure for the lower straps.
5. Adjust the straps so the leg rest sits **flush** between the footrest bars. This will allow more room on the footrests for the resident's feet.

NOTE: Extenders are available in 1", 2", 3", and 4" thicknesses and may be purchased separately. Call Customer Service at 800-447-6739.



6351

**J.T. Posey Company** 5635 Peck Road • Arcadia, CA 91006-0020 USA • Tel: 800-447-6739 or 626-443-3143 • Fax: 800-767-3933 or 626-443-5014

EU AUTH. PERSON: (MDD 93/42/EEC) MDSS GmbH, Burckhardtstr. 1, D-30163, Hannover, Germany

© 1996 J.T. Posey Company. All rights reserved.



M5015 071598