



Posey® Deluxe Guard Rail Pads

Quality Products since 1937

BED SAFETY

These practical safety aids offer maximum rail cushioning and coverage. Also available with a full-length window (#5605W, 5605WH) to help caregivers continuously monitor patients, while allowing patients to see out and feel less confined. Made with durable, firm foam; covered with soft, soil-resistant vinyl. Hook and loop straps allow fast and easy application. Rail pads wipe clean with soap & water, or disinfectant. Meets CA #117 flame retardancy standards. One pair per package.

- REF 5605** 72"L x 14"H **REF 5605W** 72"L x 14"H **REF 5606** 60"L x 14"H
- REF 5605H** 35"L x 14"H **REF 5605WH** 35"L x 14"H **REF 5607** 43"L x 14"H


APPLICATION INSTRUCTIONS:

1. Raise the side rails to the "up" position. Separate the hook and loop straps.
2. Position the side rail pad on the inside of the side rail with the straps facing out away from the bed. The top strap section should hang over the side rail.
3. Pull each strap tight around the side rail and secure back onto itself.
4. Wrap the quick-release buckles around the vertical side rail bar and secure (Model 5605W and 5605WH only).

⚠ WARNING The Deluxe Guard Rail Pads will not support weight nor prevent patients who attempt to climb over the side rail from side rail entanglement. Consult your facility restraint policy regarding patients who attempt to climb over side rails.

CLEANING INSTRUCTIONS:

Wipes Clean with mild liquid disinfectant.



J.T. Posey Company • 5635 Peck Road, Arcadia, CA 91006-0020 USA
Phone: 800-447-6739 • Fax: 800-767-3933 • Web: www.posey.com



MDSS GmbH
Schiffgraben 41
D-30175 Hannover, Germany
© 2007 J.T. Posey Company. All rights reserved.



M1302 040307



Posey® Deluxe Guard Rail Pads

Quality Products since 1937

BED SAFETY

These practical safety aids offer maximum rail cushioning and coverage. Also available with a full-length window (#5605W, 5605WH) to help caregivers continuously monitor patients, while allowing patients to see out and feel less confined. Made with durable, firm foam; covered with soft, soil-resistant vinyl. Hook and loop straps allow fast and easy application. Rail pads wipe clean with soap & water, or disinfectant. Meets CA #117 flame retardancy standards. One pair per package.

- REF 5605** 72"L x 14"H **REF 5605W** 72"L x 14"H **REF 5606** 60"L x 14"H
- REF 5605H** 35"L x 14"H **REF 5605WH** 35"L x 14"H **REF 5607** 43"L x 14"H

APPLICATION INSTRUCTIONS:

1. Raise the side rails to the "up" position. Separate the hook and loop straps.
2. Position the side rail pad on the inside of the side rail with the straps facing out away from the bed. The top strap section should hang over the side rail.
3. Pull each strap tight around the side rail and secure back onto itself.
4. Wrap the quick-release buckles around the vertical side rail bar and secure (Model 5605W and 5605WH only).

⚠ WARNING The Deluxe Guard Rail Pads will not support weight nor prevent patients who attempt to climb over the side rail from side rail entanglement. Consult your facility restraint policy regarding patients who attempt to climb over side rails.

CLEANING INSTRUCTIONS:

Wipes Clean with mild liquid disinfectant.



J.T. Posey Company • 5635 Peck Road, Arcadia, CA 91006-0020 USA
Phone: 800-447-6739 • Fax: 800-767-3933 • Web: www.posey.com



MDSS GmbH
Schiffgraben 41
D-30175 Hannover, Germany
© 2007 J.T. Posey Company. All rights reserved.



M1302 040307