



# StaySafe™ Bed

Where Comfort Meets Safety

## Posey® StaySafe™ Bed Canopy Cover/Skirt

Designed for use with the Posey StaySafe Bed to create a more child-friendly appearance.

- REF 8003P** Pediatric Cover for 8040, 8060 models
- REF 8003PS** Pediatric Skirt for 8040, 8060 models

### Application Instructions - Cover:

1. Drape the canopy cover over both 'A' frames and top center piece of the Posey StaySafe Bed (Fig. 1).
2. Make sure that the top corners are covered as shown (Fig. 2).
3. Attach the hook-and-loop tabs on each of the four corners of the canopy (Fig. 3).

**NOTE:** Remove cover before laundering the Posey StaySafe Bed canopy.

### Application Instructions - Skirt:

1. The Posey StaySafe Canopy Skirt comes in two pieces. Start at either the head or foot of the bed by snapping the two pieces together in place. You may use either snap depending on the level of tightness you want around the bottom of the StaySafe Bed Canopy Skirt (Fig. 4).
2. Drape each piece along the bottom of the StaySafe Bed Canopy Skirt and snap together at the head and foot of the bed. **NOTE:** Remove the Skirt from StaySafe Bed Canopy Skirt before laundering.

### Laundering Instructions:



### Storage And Handling

This device is designed for use in normal indoor environments. This device may be stored in ambient warehouse temperatures at normal humidity levels. Avoid excess moisture or high humidity that may damage product materials.




Fig. 1



Fig. 2



Fig. 3




Fig. 4

### Disposal

**WARNING** Properly dispose of the product per facility's policy for BIOHAZARDOUS materials.

Posey Company • 5635 Peck Road, Arcadia, CA 91006-0020 USA • www.posey.com  
Phone: 1.800.447.6739 or 1.626.443.3143 • Fax: 1.800.767.3933 or 1.626.443.5014

© 2010 Posey Company. All rights reserved.

M6182 REV B 042810



# StaySafe™ Bed

Where Comfort Meets Safety

## Posey® StaySafe™ Bed Canopy Cover/Skirt

Designed for use with the Posey StaySafe Bed to create a more child-friendly appearance.

- REF 8003P** Pediatric Cover for 8040, 8060 models
- REF 8003PS** Pediatric Skirt for 8040, 8060 models

### Application Instructions - Cover:

1. Drape the canopy cover over both 'A' frames and top center piece of the Posey StaySafe Bed (Fig. 1).
2. Make sure that the top corners are covered as shown (Fig. 2).
3. Attach the hook-and-loop tabs on each of the four corners of the canopy (Fig. 3).

**NOTE:** Remove cover before laundering the Posey StaySafe Bed canopy.

### Application Instructions - Skirt:

1. The Posey StaySafe Canopy Skirt comes in two pieces. Start at either the head or foot of the bed by snapping the two pieces together in place. You may use either snap depending on the level of tightness you want around the bottom of the StaySafe Bed Canopy Skirt (Fig. 4).
2. Drape each piece along the bottom of the StaySafe Bed Canopy Skirt and snap together at the head and foot of the bed. **NOTE:** Remove the Skirt from StaySafe Bed Canopy Skirt before laundering.

### Laundering Instructions:



### Storage And Handling

This device is designed for use in normal indoor environments. This device may be stored in ambient warehouse temperatures at normal humidity levels. Avoid excess moisture or high humidity that may damage product materials.




Fig. 1



Fig. 2




Fig. 3




Fig. 4

### Disposal

**WARNING** Properly dispose of the product per facility's policy for BIOHAZARDOUS materials.

Posey Company • 5635 Peck Road, Arcadia, CA 91006-0020 USA • www.posey.com  
Phone: 1.800.447.6739 or 1.626.443.3143 • Fax: 1.800.767.3933 or 1.626.443.5014

© 2010 Posey Company. All rights reserved.

M6182 REV B 042810