



Posey® Elbow Splints

NICU/PEDS

- An external hand, wrist, elbow and knee splint with a quilted cover and inner lining.
- Adjustable hook closure. Elbow Splints may be pinned to the child's clothing to hold them in place. One per package.

REF 4642 - Infant (2 years), 4"L (10 cm)

REF 4643 - Pediatric (2-7 years), 6¾"L (17 cm)

APPLICATION INSTRUCTIONS

1. Release the hook-and-loop adjustment straps.
2. Wrap the Elbow Splints around the arm and center over the elbow.
3. Secure the hook fastener to the quilted material. Leave enough room to easily insert two fingers between the device and the patient's limb to maintain adequate circulation.

WARNING Monitor skin condition frequently. DO NOT overtighten or impair circulation!

NEVER alter or repair this product. ALWAYS inspect before each use: Check for broken stitches or parts; torn, cut or frayed material; or hook and loop fasteners that do not hold securely. DO NOT use soiled or damaged products. Doing so may result in serious injury or death. Dispose of damaged products per facility policy for BIOHAZARDOUS material.

WARNING DO NOT allow patients to ingest product material.

STORAGE AND HANDLING

This device is designed for use in normal indoor environments. This device may be stored in ambient warehouse temperatures at normal humidity levels. Avoid excess moisture or high humidity that may damage product materials.



LAUNDERING INSTRUCTIONS:



Remove tongue depressor stiffeners before laundering.



Posey Company • 5635 Peck Road, Arcadia, CA 91006-0020 USA • www.posey.com
Phone: 1.800.447.6739 or 1.626.443.3143 • Fax: 1.800.767.3933 or 1.626.443.5014



MDSS GmbH
Schiffgraben 41
D-30175 Hannover, Germany



© 2009 Posey Company. All rights reserved.

M6230 072309



Posey® Elbow Splints

NICU/PEDS

- An external hand, wrist, elbow and knee splint with a quilted cover and inner lining.
- Adjustable hook closure. Elbow Splints may be pinned to the child's clothing to hold them in place. One per package.

REF 4642 - Infant (2 years), 4"L (10 cm)

REF 4643 - Pediatric (2-7 years), 6¾"L (17 cm)

APPLICATION INSTRUCTIONS

1. Release the hook-and-loop adjustment straps.
2. Wrap the Elbow Splints around the arm and center over the elbow.
3. Secure the hook fastener to the quilted material. Leave enough room to easily insert two fingers between the device and the patient's limb to maintain adequate circulation.

WARNING Monitor skin condition frequently. DO NOT overtighten or impair circulation!

NEVER alter or repair this product. ALWAYS inspect before each use: Check for broken stitches or parts; torn, cut or frayed material; or hook and loop fasteners that do not hold securely. DO NOT use soiled or damaged products. Doing so may result in serious injury or death. Dispose of damaged products per facility policy for BIOHAZARDOUS material.

WARNING DO NOT allow patients to ingest product material.

STORAGE AND HANDLING

This device is designed for use in normal indoor environments. This device may be stored in ambient warehouse temperatures at normal humidity levels. Avoid excess moisture or high humidity that may damage product materials.



LAUNDERING INSTRUCTIONS:



Remove tongue depressor stiffeners before laundering.



Posey Company • 5635 Peck Road, Arcadia, CA 91006-0020 USA • www.posey.com
Phone: 1.800.447.6739 or 1.626.443.3143 • Fax: 1.800.767.3933 or 1.626.443.5014



MDSS GmbH
Schiffgraben 41
D-30175 Hannover, Germany



© 2009 Posey Company. All rights reserved.

M6230 072309