

Are You Protecting Your Patients from Falls?



6020, 6021, 6023, 6024
6025, 6026, 6027, 6027R


Posey Floor Cushions help reduce the risk of fall-related injuries by absorbing up to 85%* of the impact from a fall. Available in standard and Beveled-Edge models, the low-profile Floor Cushions feature wipe-clean vinyl covers with non-skid bottoms for easier ambulation.



6020



6023



6027



6026

We have more practical ways to protect your patients.

Models 6020, 6021, 6023, 6025, 6026, 6027, 6027R

- EVA foam absorbs up to 85% of impact forces from a fall.
- One-inch thick foam creates low profile edge.
- Non-skid bottom.
- Tri-fold design with convenient carrying handle.
- Meets CA #117, Section E Flame Retardancy Standards.

Model 6024

- Helps absorb up to 65% of the impact from a fall.
- High density, polyurethane foam.
- Folds in half with convenient carrying handle.

For additional Fall Protection, Posey Floor Cushions may be used with the 8250 Exit Alarm Mat and any sensor-compatible Posey Fall-Management Monitors.

NOTE: Due to the random possibility of falls, the Posey Company makes no guarantee, express or implied, that the user is protected from hip trauma.

*All models tested according to ASTM test protocols for shock absorbency.

*Using a Posey Beveled Floor Cushion decreased the risk of injury to less than 1% for bedrail height for feet-first falls.**

* Bowers R, Lloyd J, et. al. Biomechanical Evaluation of Injury Severity Associated with Patient Falls from Bed. *Rehabilitation Nursing* 2008, Vol. 33, No. 6.

To request a Posey Catalog or speak with an associate, call 1.800.44.POSEY (76739) or visit us online at Posey.com

Posey Company • 5635 Peck Road, Arcadia, CA 91006-0020 USA • Phone: 1.800.447.6739 or 1.626.443.3143
Fax: 1.800.767.3933 or 1.626.443.5014 • www.posey.com ©2010 Posey Company. All rights reserved. MK1245 REV C 101610





Posey Beveled Floor Cushions

- Features beveled edges on three sides, which helps reduce tripping hazards and accommodates walkers and wheelchairs.
- Narrow Floor Cushion is beveled on all four sides.
- Made of low profile EVA foam that absorbs up to 85% of impact forces from a fall.
- Tested according to ASTM test protocols for shock absorbency. Covered in wipe-clean vinyl and folds into thirds for easy storage with hook-and-loop straps.
- Posey Glow model includes a luminescent border that glows for a limited time after exposure to light.

DIMENSIONS:

6023 Beveled Floor Cushion	Flat: 70”L x 38”W x 1”H, Folded: 24”L x 38”W x 3”H
6025 Beveled Floor Cushion, Wide	Flat: 70”L x 44”W x 1”H, Folded: 24”L x 44”W x 3”H
6026 Beveled Floor Cushion, Narrow	Flat: 70”L x 29”W x 1”H, Folded: 24”L x 29”W x 3”H
6027 Beveled Floor Cushion, Premium	Flat: 70”L x 38”W x 1”H, Folded: 24”L x 38”W x 3”H
6027R Beveled Floor Cushion, Posey Glow	Flat: 70”L x 38”W x 1”H, Folded: 24”L x 38”W x 3”H



Posey Floor Cushions

- EVA foam absorbs up to 85% of impact forces from a fall
- Tested according to ASTM test protocols for shock absorbency. Non-skid bottom helps prevent slipping and aids in patient ambulation.
- Bacteriostatic vinyl covering is available in either blue or mauve, for easier visibility in dimly lit rooms.
- 1” thick foam creates a low profile edge that is easy for patients to step on and off.
- Provides full length coverage of the bed exit area and easily slides under the bed when not in use. Helps reduce injury from falls when side rails are down.
- Tri-fold design and convenient carrying handle allow easy portability and storage.
- Wipes clean with liquid disinfectant.
- Meets FMVSS302 Flame Retardancy Standard.

DIMENSIONS:

6020 Floor Cushion	Flat: 67”L x 26”W x 1”H, Folded: 23”L x 26”W x 3”H
6021 Two-Tone Floor Cushion	Flat: 67”L x 26”W x 1”H Folded: 23”L x 26”W x 3”H



Posey Economy Floor Cushion

- Made of high-density polyurethane, impact-absorbing foam.
- Helps reduce the possibility of impact related injuries when side rails are down.
- Provides full coverage of the bed exit area and conveniently folds in half.
- Low profile pad absorbs up to 65% of impact forces from a fall and the vinyl cover wipes clean with liquid disinfectant.
- Foam meets CA #117 Section E Flame Retardancy Standards, while the vinyl cover meets California Title 19 Flame Retardancy Standards.

6024 Economy Floor Cushion

DIMENSIONS:

Flat: 72”L x 24”W x 1½”H, Folded: 36”L x 24”W x 3”H

To request a Posey Catalog or speak with an associate, call 1.800.44.POSEY (76739) or visit us online at Posey.com

Posey Company • 5635 Peck Road, Arcadia, CA 91006-0020 USA • Phone: 1.800.447.6739 or 1.626.443.3143
Fax: 1.800.767.3933 or 1.626.443.5014 • www.posey.com

