



# Posey®



# PRO-heelLx®

Prevention on. Pressure off.



The outer fabric is made of breathable nylon, allowing moisture and heat to migrate away from the foot.

**The Posey PRO-heelLx®** is used for offloading the heel, and can be used on either the left or right foot.



2" heel lift helps isolate the heel.

### WEDGE - 6219

Optional Wedge helps prevent foot rotation

The PRO-heelLx® comes in 3 sizes to customize fit any patient type:

**LARGE - 6218L**  
(for sizes over men's 11)

**MEDIUM - 6218**  
(fits up to men's size 11)

**SMALL - 6218S**  
(fits women's shoe size up to 6)



Size	PRODUCT DIMENSIONS			SUGGESTED FITTING GUIDELINES		
	Horizontal Calf Section	Vertical Foot Section	Side Panels	Calf Circumference	Shoe Sizes	
					Women	Men
<b>6218S</b> - Small	4½" W x 11" L (11 cm x 28 cm)	4½" W x 9" H (11 cm x 23 cm)	4" H (10 cm)	Up to 12" (30 cm)	Up to 6	-
<b>6218</b> - Medium	5½" W x 13" L (14 cm x 33 cm)	5½" W x 10" H (14 cm x 25 cm)	5" H (13 cm)	Up to 16" (41 cm)	7-11	Up to 11
<b>6218L</b> - Large	6½" W x 16" L (17 cm x 41 cm)	6½" W x 13" H (17 cm x 33 cm)	5" H (13 cm)	Up to 20" (51 cm)	-	12+

To request a Posey Catalog or speak with an associate, call 1.800.44.POSEY (76739) or visit us online at Posey.com

Posey Company • 5635 Peck Road, Arcadia, CA 91006-0020 USA • Phone: 1.800.447.6739 or 1.626.443.3143  
Fax: 1.800.767.3933 or 1.626.443.5014 • www.posey.com ©2012 Posey Company. All rights reserved. SA0411A REV A 082912

