



# Posey Company

5635 Peck Road

Arcadia, CA 91006-0020 USA

Toll Free: (800) 44-POSEY or (800) 447-6739

Web: [www.posey.com](http://www.posey.com)

**FAQ No.:** 09

**FAQ:**

What temperatures are associated with our symbols "Dry Low" and "Dry High".

---

The Posey Company recommends the following standards set by the National Cleaners Association New York, New York.

"Dry Low" means a dryer temperature of 100 degrees F to 105 degrees F.

"Dry High" means a dryer temperature of 140 degrees F to 180 degrees F.

For complete cleaning and use instructions for all Posey products, please refer to package insert instructions accompanying each product. Copies of all package insert instructions are available on the Posey website at [www.posey.com](http://www.posey.com).

We hope that this answers your questions. If you have any further questions, please feel free to call Posey Customer Service at (800) 44-Posey or (800) 447-6739.