



Posey Company

5635 Peck Road

Arcadia, CA 91006-0020 USA

Toll Free: (800) 44-POSEY or (800) 447-6739

Web: www.posey.com

FAQ No.: 31

FAQ:

Can the Posey Belt #4125, Posey Padded Belts #4127 and #4135 be used on patients in chairs other than a wheelchair ?

Yes, if the belt can be applied according to the instructions accompanying each product. These instructions require the belt to be applied at a 45 degree angle over the hips, and secured at a juncture of the frame. This helps prevent the patient from sliding under the belt and suspending themselves off the floor. If the belt is not attached properly, and the patient is able to slide off the seat and suspend themselves off the floor, they may be unable to call for help and asphyxiate.

For complete application and care instructions for all Posey products, please refer to package insert accompanying each Posey product.

We hope this answers your questions. If you have any further questions, please feel free to call Posey Customer Service at (800) 44-POSEY or (800) 447-6739.